



SonoSite Transducer Technology

The **MicroMaxx® system** offers a wide range of transducers for multiple clinical applications. These transducers are designed with broadband signal processing, which provides access to more information by transmitting and receiving through a broader bandwidth of frequencies.* This technology allows a single transducer to image over a greater range of depths, expanding clinical utility. Unique in the industry, SonoSite transducers yield an uptime of greater than 99.9% per year.



L38e/10-5

38-mm broadband linear array

Major applications

Small parts, breast, vascular, nerve, IMT, musculoskeletal, superficial
Biopsy guide available

HFL38/13-6

38-mm broadband linear array

Major applications

Breast, small parts, nerve, vascular, IMT, musculoskeletal
Biopsy guide available

L25e/13-6

25-mm broadband linear array

Major applications

Nerve, musculoskeletal, vascular and superficial
Transverse biopsy guide available

SLA/13-6

26-mm broadband linear array

Major applications

Vascular, musculoskeletal, superficial, nerve

C60e/5-2

60-mm broadband curved array

Major applications

Abdominal, obstetrics, gynecology
Biopsy guide available

C11e/8-5

11-mm broadband curved array

Major applications

Abdominal, nerve



ICT/8-5

11-mm broadband tightly curved array

Major applications

Obstetrics, gynecology
Biopsy guide available



P17/5-1

17-mm broadband phased array

Major applications

Cardiac, abdominal, obstetrics
Biopsy guide available



P10/8-4

10-mm broadband phased array

Major applications

Pediatric cardiology, pediatric abdominal, neonatal heads, nerve
Biopsy guide available



D2/2

Pedof CW Pencil Transducer

Major application

Cardiac



TEE/8-3

Multiplane transesophageal
180-degree rotation of the imaging plane, providing a 360-degree field of view

Major applications

Adult transesophageal imaging with full cardiac calculation package

*With the exception of the D2/2, a single frequency transducer.
©2006 SonoSite, Inc. All rights reserved. Subject to change. MKT00767 Rev C 9-06

WORLDWIDE HEADQUARTERS

SonoSite, Inc.
21919 30th Drive SE, Bothell, WA 98021-3904 USA
Telephone: (425) 951-1200 Toll Free: (877) 657-8050 Fax: (425) 951-6800
E-mail: admin@sonosite.com Web: www.sonosite.com